

INCIDENT INFORMATION			
PATIENT'S NAME	CHECKLIST COMPLETED BY (PROVIDER)	INCIDENT DATE (DD / MM / YY)	TIME (24 hr)

When the listed Indications are met, Contraindications ruled out and Conditions achieved, a Qualified Provider¹ may apply an AutoPulse as an adjunct to manual CPR on adult non-traumatic cardiac arrest patients who require ongoing resuscitation.

INDICATIONS FOR USE

- Adult (>18 years of age) patients with medical (non-traumatic) cardiac arrest

CONTRAINDICATIONS OF USE

1. Patients sustaining cardiac arrest due to trauma
2. Pediatric patients (<18 years of age)
3. Patient's physical parameters that are not within the AutoPulse's operating parameters:

AUTOPULSE OPERATING PARAMETERS	
PATIENT PARAMETERS	AUTOPULSE SPECIFICATIONS
Patient chest circumference permitted	29.9" to 51.2" / 76 cm to 130 cm
Patient chest width permitted	9.8" to 15" / 25 cm to 38 cm
Maximum patient weight permitted	300 lbs. / 136 kg

CONDITIONS OF USE (ALL CONDITIONS MUST BE ACHIEVED PRIOR TO APPLICATION OF AN AUTOPULSE)

- Completed a Primary Assessment
- Provided and maintained all necessary critical interventions
- Activated the local emergency response system and notified EMS
- Conducted a targeted Past Medical History if possible
- Ruled out contraindications
- Patient must be observed by a Qualified Responder or higher medically qualified person during transport
- Application of an AutoPulse is to be applied as per PEAK's APTP Training Manual and AutoPulse Operational Guideline (separate document)

¹ Qualified Provider - an individual who has received specific PEAK training.