



# *EMS Vacuum Immobilisation Catalogue 2023*



# **RedVac EMS Vacuum Immobilisation**

## **Content Page**

**RedVac made in Austria 3**

**The RedVac System 4-5**

**Mattresses 6-10**

**Head Immobilisation 11**

**Extrication Device 11**

**Extremity Splints 12-18**

**Pediatric 19-20**

**Pumps, Bags & Accessories 21**

**Thermo & HRS Rescue Sacks 22**

**Helicopter Rescue Sack 23**



# Developed and manufactured in Austria

Kohlbrat & Bunz products are developed and manufactured 900 meters above sea level in the middle of the Austrian Alps in its own facilities by its experienced and meticulously working staff.


People living in the Alps know very well that products used under the extreme conditions of alpine rescue have to work perfectly.

To meet these rigorous demands our products are developed in close cooperation with the mountain rescue teams and other demanding emergency services.

This cooperation has brought about a selection of the best materials and optimal product design leading to products of outstanding performance and quality. To insure the best possible quality Kohlbrat & Bunz work in accordance with EN ISO 9001 & EN ISO 13485.

## Worldwide in use

We are very proud of the confidence shown to us and our products by many well-known and demanding emergency services worldwide and are happy to show some examples of their beautiful logos that they have asked us to put on the RedVac products we have made for them.





# The RedVac System

for **PATIENT SAFETY, PERFECT HYGIENE & EASE OF USE**

- **The RedVac CLICK-System**  
for attachment of straps and handles

*The unique CLICK-System makes removal of straps and handles for repositioning, replacement, cleaning and disinfection very easy!*



- **The RedVac Multi-Chamber-Systems**  
- provides unmatched **EASE OF USE**

*The Multi-Chamber-Systems keeps the internal beads evenly distributed eliminating the need for difficult and time consuming manual spreading of the beads. This enables fast, safe and easy application in all situations.*

**The RedVac Multi-Chamber-System!**

*The internal sack system keeps the beads evenly distributed in up to 120 chambers.*



**The RedVac DOT-Chamber-System!**

*The layers of the outer cover are welded together dot by dot creating chambers between the dots keeping the internal beads in their correct position, evenly distributed.*



- **Optional attachment systems for straps and handles**

**Cleats**



**Hook & Loop**



**Red Stripes**



- **Robust & durable exterior Material**

To achieve outstanding durability, abrasion resistance, tightness, easy cleaning and disinfection, comfort and safety for the patient special materials were chosen. It makes it possible to use the RedVac products in a wide temperature range, from -40°C up to +80°C.



- **Restraint Straps, made of polypropylene (PP) or wipe-clean coated fabric**

Colour coded patient restraint straps, adjustable from both sides, easy to clean, reposition or replace.



- **Vacuum Valve, screw-on, replaceable**

The valve is robust and easy to use, even when wearing gloves. It closes automatically when suction stops preventing loss of vacuum if slipping with the pump. Comes with suitable adapter for most vacuum pumps.



- **Vacuum Valve, welded-on, non-replaceable**

The valve is robust and easy to use, even when wearing gloves. It closes automatically when suction stops preventing loss of vacuum if slipping with the pump.



- **Label System**

An indestructible label with individual serial number and date of manufacture is welded onto every product. Optional labels with property indicator and distributor logo are available.



- **X-ray, CT and MRI compatible**



- **CE**

All RedVac products are CE marked and conform with the European directive 93/42/EEC for medical products class 1.




- **ISO9001 - ISO13485 certified by**




# RedVac EMS Vacuum-Mattresses




## Body shape

	<b>VM01X02</b> <b>VM01HX02</b>	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	<b>VM01CX02</b> <b>VM01CHX02</b>	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 40-88-50cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 4		


## Body shape

	<b>VM03X03</b> <b>VM03HX03</b>	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	<b>VM03CX03</b> <b>VM03CHX03</b>	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 50-88-50cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 4		



## Body shape. Slits for gurney safety belts in neck area.

	<b>VM03X05</b> <b>VM03HX05</b>	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	<b>VM03CX05</b> <b>VM03CHX05</b>	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 50-98-50cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 4		

## Body shape. Slits for gurney safety belts in neck area. Glass fiber rod for cervical protection.



	<b>VM03X07</b> <b>VM03HX07</b>	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	<b>VM03CX07</b> <b>VM03CHX07</b>	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 50-88-50cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 4		

## Body shape - compact size for helicopter transportation


	<b>VM08X06</b> <b>VM08HX06</b>	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	<b>VM08CX06</b> <b>VM08CHX06</b>	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 40-78-40cm Length: 195cm Weight: 6kg Carry handles: 8 Patient restraint straps: 4		
	<b>VM09X01</b> <b>VM09HX01</b>	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	<b>VM09CX01</b> <b>VM09CHX01</b>	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 40-60-40cm Length: 195cm Weight: 5kg Carry handles: 8 Patient restraint straps: 4		

## Body shape with head fixation.


### Slits for gurney safety belts in neck area. Glass fiber rod for cervical protection.

	<b>VM200X06</b> <b>VM200HX06</b>	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	<b>VM200CX06</b> <b>VM200CHX06</b>	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 40-88-60cm Length: 200cm Weight: 7kg Carry handles: 8 Patient restraint straps: 6		
	<b>VM201X10</b> <b>VM201X10</b>	CLEATS with PP restraint straps CLEATS with wipe-clean restraint straps
	<b>VM201CX10</b> <b>VM201CHX10</b>	CLICK with PP restraint straps CLICK with wipe-clean restraint straps
Width head/mid/foot: 40-100-60cm Length: 210cm Weight: 8kg Carry handles: 8 Patient restraint straps: 6		


**Head fixation and slits in neck area for gurney safety belts.**

	<p><b>VM500X03</b> CLEATS with PP restraint straps  <b>VM500HX03</b> CLEATS with wipe-clean restraint straps  <b>VM500CX03</b> CLICK with PP restraint straps  <b>VM500CHX03</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 60-88-88cm Length: 200cm Weight: 7,5kg            Carry handles: 8 Patient restraint straps: 5</p>
	<p><b>VM501X03</b> CLEATS with PP restraint straps  <b>VM501HX03</b> CLEATS with wipe-clean restraint straps  <b>VM501CX03</b> CLICK with PP restraint straps  <b>VM501CHX03</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 60-100-100cm Length: 210cm Weight: 8,5kg            Carry handles: 12 Patient restraint straps: 5</p>
	<p><b>VM502X01</b> CLEATS with PP restraint straps  <b>VM502HX01</b> CLEATS with wipe-clean restraint straps  <b>VM502CX01</b> CLICK with PP restraint straps  <b>VM502CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 60-80-80cm Length: 200cm Weight: 7kg            Carry handles: 8 Patient restraint straps: 5</p>

**Head fixation and slits in neck area for gurney safety belts. Glassfiber rod for cervical protection.**



	<p><b>VM5000X01</b> CLEATS with PP restraint straps  <b>VM5000HX01</b> CLEATS with wipe-clean restraint straps  <b>VM5000CX01</b> CLICK with PP restraint straps  <b>VM5000CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 65-88-88cm Length: 200cm Weight: 7kg            Carry handles: 8 Patient restraint straps: 5</p>
	<p><b>VM5001X01</b> CLEATS with PP restraint straps  <b>VM5001HX01</b> CLEATS with wipe-clean restraint straps  <b>VM5001CX01</b> CLICK with PP restraint straps  <b>VM5001CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 65-100-100cm Length: 210cm Weight: 8kg            Carry handles: 12 Patient restraint straps: 5</p>
	<p><b>VM5002X01</b> CLEATS with PP restraint straps  <b>VM5002HX01</b> CLEATS with wipe-clean restraint straps  <b>VM5002CX01</b> CLICK with PP restraint straps  <b>VM5002CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 65-80-80cm Length: 200cm Weight: 6kg            Carry handles: 8 Patient restraint straps: 5</p>

**Body shape with head fixation. Slits for gurney safety belts. Glass fiber rod for cervical protection. Additional outer protective layer.**


	<p><b>VM7000X03</b> CLEATS with PP restraint straps  <b>VM7000HX03</b> CLEATS with wipe-clean restraint straps  <b>VM7000CX03</b> CLICK with PP restraint straps  <b>VM7000CHX03</b> CLICK with wipe-clean restraint straps</p> <p>Width head-middle-foot: 60-130-60cm Length: 205cm Weight: 10kg            Carry handles, integrated in protective layer: 10 Carry handles removable: 4            Patient restraint straps: 5</p>
	<p><b>VM7001CX04</b> CLICK with PP restraint straps  <b>VM7001CHX04</b> CLICK with wipe-clean restraint straps</p> <p>Width head-middle-foot: 60-110-60cm Length: 205cm Weight: 9,5kg            Carry handles, integrated in protective layer: 10            Patient restraint straps: 5</p>
	<p><b>VM7002CFX01</b> CLICK with PP restraint straps  <b>VM7002CFHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width head-middle-foot: 60-110-60cm Length: 205cm Weight: 9,5kg            Carry handles removable: 14 Patient restraint straps: 7</p>



**Head fixation. Slits for gurney safety belts in neck area. Glassfiber rod for cervical protection**



	<p><b>VM11X04</b> CLEATS with PP restraint straps  <b>VM11HX04</b> CLEATS with wipe-clean restraint straps  <b>VM11CX04</b> CLICK with PP restraint straps  <b>VM11CHX04</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 80-88-80cm Length: 200cm Weight: 7kg            Carry handles: 8 Restraint straps: 4</p>
	<p><b>VM1101X01</b> CLEATS with PP restraint straps  <b>VM1101HX01</b> CLEATS with wipe-clean restraint straps  <b>VM1101CX01</b> CLICK with PP restraint straps  <b>VM1101CHX04</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 80-100-80cm Length: 210cm Weight: 8kg            Carry handles: 12 Restraint straps: 4</p>
	<p><b>VM1111X01</b> CLEATS with PP restraint straps  <b>VM1111HX01</b> CLEATS with wipe-clean restraint straps  <b>VM1111CX01</b> CLICK with PP restraint straps  <b>VM1111CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 80-108-80cm Length: 210cm Weight: 8,5kg            Carry handles: 12 Restraint straps: 4</p>
	<p><b>VM1112X01</b> CLEATS with PP restraint straps  <b>VM1112HX01</b> CLEATS with wipe-clean restraint straps  <b>VM1112CX01</b> CLICK with PP restraint straps  <b>VM1112CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 80-108-80cm Length: 220cm Weight: 8,5kg            Carry handles: 8 Restraint straps: 4</p>

**Rectangular shape. Slits for gurney safety belts in the neck area Glassfiber rod for cervical protection.**


	<p><b>VM11620X02</b> CLEATS with PP restraint straps  <b>VM11620HX02</b> CLEATS with wipe-clean restraint straps  <b>VM11620CX02</b> CLICK with PP restraint straps  <b>VM11620CHX02</b> CLICK with wipe-clean restraint straps</p> <p>Width: 60cm Length: 200cm Weight: 5kg            Carry handles: 8 Restraint straps: 4</p>
	<p><b>VM11720X01</b> CLEATS with PP restraint straps  <b>VM11720X01</b> CLEATS with wipe-clean restraint straps  <b>VM11720CX01</b> CLICK with PP restraint straps  <b>VM11720CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width: 70cm Length: 200cm Weight: 6kg            Carry handles: 8 Restraint straps: 4</p>
	<p><b>VM11820X01</b> CLEATS with PP restraint straps  <b>VM11820X01</b> CLEATS with wipe-clean restraint straps  <b>VM11820CX01</b> CLICK with PP restraint straps  <b>VM11820CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width: 80cm Length: 200cm Weight: 7kg            Carry handles: 8 Restraint straps: 4</p>
	<p><b>VM11920X01</b> CLEATS with PP restraint straps  <b>VM11920HX01</b> CLEATS with wipe-clean restraint straps  <b>VM11920CX01</b> CLICK with PP restraint straps  <b>VM11920CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width: 90cm Length: 200cm Weight: 7,5kg            Carry handles: 8 Restraint straps: 4</p>
	<p><b>VM1201X01</b> CLEATS with PP restraint straps  <b>VM1201HX01</b> CLEATS with wipe-clean restraint straps  <b>VM1201CX01</b> CLICK with PP restraint straps  <b>VM1201CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width: 100cm Length: 210cm Weight: 8,5kg            Carry handles: 12 Restraint straps: 4</p>



**Tapering shape at head and foot area**


	<p><b>VM10X02</b> CLEATS with PP restraint straps  <b>VM10HX02</b> CLEATS with wipe-clean restraint straps  <b>VM10CX02</b> CLICK with PP restraint straps  <b>VM10CHX02</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 80-88-80cm Length: 200cm Weight: 7kg            Carry handles: 8 Restraint straps: 4</p>
	<p><b>VM8021X01</b> CLEATS with PP restraint straps  <b>VM8021HX01</b> CLEATS with wipe-clean restraint straps  <b>VM8021CX01</b> CLICK with PP restraint straps  <b>VM8021CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 80-100-80cm Length: 210cm Weight: 8,5kg            Carry handles: 12 Restraint straps: 4</p>
	<p><b>VM8121X01</b> CLEATS with PP restraint straps  <b>VM8121HX01</b> CLEATS with wipe-clean restraint straps  <b>VM8121CX01</b> CLICK with PP restraint straps  <b>VM8121CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width head/mid/foot: 80-108-80cm Length: 220cm Weight: 9kg            Carry handles: 12 Restraint straps: 4</p>

**Rectangular shape**

	<p><b>VM620X01</b> CLEATS with PP restraint straps  <b>VM620HX01</b> CLEATS with wipe-clean restraint straps  <b>VM620CX01</b> CLICK with PP restraint straps  <b>VM620CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width: 60cm Length: 200cm Weight: 6kg            Carry handles: 8 Restraint straps: 4</p>
	<p><b>VM720X01</b> CLEATS with PP restraint straps  <b>VM720HX01</b> CLEATS with wipe-clean restraint straps  <b>VM720CX01</b> CLICK with PP restraint straps  <b>VM720CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width: 70cm Length: 200cm Weight: 6,5kg            Carry handles: 8 Restraint straps: 4</p>
	<p><b>VM820X01</b> CLEATS with PP restraint straps  <b>VM820HX01</b> CLEATS with wipe-clean restraint straps  <b>VM820CX01</b> CLICK with PP restraint straps  <b>VM820CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width: 80cm Length: 200cm Weight: 7kg            Carry handles: 8 Restraint straps: 4</p>
	<p><b>VM1021X01</b> CLEATS with PP restraint straps  <b>VM1021HX01</b> CLEATS with wipe-clean restraint straps  <b>VM1021CX01</b> CLICK with PP restraint straps  <b>VM1021CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width: 100cm Length: 210cm Weight: 8,5kg            Carry handles: 12 Restraint straps: 4</p>
<p><b>VM1121X01</b> CLEATS with PP restraint straps  <b>VM1121HX01</b> CLEATS with wipe-clean restraint straps  <b>VM1121CX01</b> CLICK with PP restraint straps  <b>VM1121CHX01</b> CLICK with wipe-clean restraint straps</p> <p>Width: 108cm Length: 210cm Weight: 9kg Handles: 12            Carry handles: 12 Restraint straps: 4</p>	

**For obese patients**


**Slits for gurney's safety belts in the neck area. Glass fiber rod for cervical protection**

	<p><b>VM1421X01</b> CLEATS with PP restraint straps</p>
	<p><b>VM1421CX01</b> CLICK with PP restraint straps  <b>VM1421CHX01</b> CLICK with wipe-clean PVC restraint straps</p> <p>Width head-mid-foot: 80-140-140cm Length: 210cm Weight: 11kg          Carry handles: 10 Restraint straps: 5</p>

**Additional protective layer with integrated handles**

	<p><b>VM2001CX01</b> CLICK with PP restraint straps  <b>VM2001CHX01</b> CLICK with wipe-clean PVC restraint straps</p> <p>Width head-mid-foot: 80-100-80cm Length: 205cm Weight: 8,5kg          Carry handles: 10 Restraint straps: 4</p>
---	---

**Additional protective layer with integrated handles**



	<p><b>VM4000CX01</b> CLICK with PP restraint straps  <b>VM4000CHX01</b> CLICK with wipe-clean PVC restraint straps</p> <p>Width: 95cm Length: 200cm Weight: 8,5kg          Carry handles: 10 Restraint straps: 4</p>
---	--

**The RedVac Combi-Vacuum-Mattress**

**2 in 1 - Combines the standard mattress and the vacuum mattress**

**Body shape with head fixation. Slits for gurney safety belts.**

**Glass fiber rod for cervical protection.**

	<p><b>Combi-Vacuum-Mattress, without Zipper</b></p> <p><b>VM9001CX01</b> CLICK with PP restraint straps  <b>VM9001CHX01</b> CLICK with wipe-clean PVC restraint straps</p> <p>Width head/middle/foot: 75-104-50cm          Length: 200cm Weight: 9kg</p>
	<p><b>Combi-Vacuum-Mattress, with Zipper</b></p> <p><b>VM9002CX02</b> CLICK with PP restraint straps  <b>VM9002CHX02</b> CLICK with wipe-clean PVC restraint straps</p> <p>Width head/middle/foot: 75-104-50cm          Length: 200cm Weight: 9kg</p>
	<p><b>Vacuum Pillow</b>  <b>VK200</b></p> <p>Width: 57cm Length: 39cm Weight: 0,5kg</p>

# The RedVac Vacuum-Head-Immobiliser - for the Spine Board & the Scoop Stretcher

- Quick and easy positioning of the immobiliser
- No manipulation of head and cervical spine
- Adjusts automatically to any head shape and size
- Ears remain free for unrestricted hearing and medical examination
- Solid and comfortable immobilisation
- Vacuum valves: Removable: R-Valve, Fixed: F-Valve



## Specification

	<p><b>Head Immobiliser</b>  <b>VK98X07</b> R-Valve  <b>VK98X07V</b> F-Valve            Width: 38cm Length: 31cm            Weight: 0,7kg</p> <p><b>Base plates</b>  <b>GP08</b> for scoop stretcher  <b>GP11</b> for back board</p>
	<p><b>Head Immobiliser with buckles for attachment</b>  <b>VK94X17</b> R-Valve  <b>VK94X17V</b> F-Valve            Width: 38cm Length: 31cm            Weight: 0,7kg</p>

# RedVac Vacuum Extrication Device

- the answer to your extrication challenges!
- fits every patient, pregnant – overweight ... no problem!
- comfortable – no respiratory impairment!

## Specification

<b>VED24X02</b>	R-Valve
<b>VED24X02V</b>	F-Valve
Width, head: 40cm	
Width, middle: 106cm	
Length: 94cm	
Weight: 2,5kg	





# RedVac Vacuum Splints

Optional attachment systems for patient straps:

CLICK      HOOK & LOOP      RED STRIPES      CLEATS

Vacuumvalves:

Removable=R      Fixed=F






Patient straps: Washable of PP or Wipe-clean of PVC






## RedVac Vacuum Splints "MOZART"







MOZART splints with RED STRIPES & replaceable PP straps

	<p>Forearm Splints (Colour blue – blue)</p> <table border="0"> <tr> <td><b>VS60X01</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS60X01V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS61X01</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS61X01V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 3    Width: 35cm    Length: 50cm    Weight: 0,33kg</p>	<b>VS60X01</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber	<b>VS60X01V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber	<b>VS61X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	<b>VS61X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber
<b>VS60X01</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber											
<b>VS60X01V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber											
<b>VS61X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber											
<b>VS61X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber											
	<p>Arm Splints (Colour blue – blue)</p> <table border="0"> <tr> <td><b>VS62X01</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS62X01V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS63X01</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS63X01V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5    Width top/foot: 43-70cm    Length: 56cm    Weight: 0,9kg</p>	<b>VS62X01</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber	<b>VS62X01V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber	<b>VS63X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	<b>VS63X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber
<b>VS62X01</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber											
<b>VS62X01V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber											
<b>VS63X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber											
<b>VS63X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber											
	<p>Leg Splints (Colour blue – blue)</p> <table border="0"> <tr> <td><b>VS64X01</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS64X01V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS65X01</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS65X01V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5    Width top/foot: 64-70cm    Length: 103cm    Weight: 1,9kg</p>	<b>VS64X01</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber	<b>VS64X01V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber	<b>VS65X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	<b>VS65X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber
<b>VS64X01</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber											
<b>VS64X01V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber											
<b>VS65X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber											
<b>VS65X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber											

## MOZART splints with replaceable Loop straps

	<p>Forearm Splints (Colour blue – blue)</p> <table border="0"> <tr> <td><b>VS60X20V</b></td> <td>Hook &amp; Loop straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS61X20V</b></td> <td>Hook &amp; Loop straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 2    Width: 35cm    Length: 50cm    Weight: 0,33kg</p>	<b>VS60X20V</b>	Hook & Loop straps, F-Valve	1-Chamber	<b>VS61X20V</b>	Hook & Loop straps, F-Valve	Multi-Chamber
<b>VS60X20V</b>	Hook & Loop straps, F-Valve	1-Chamber					
<b>VS61X20V</b>	Hook & Loop straps, F-Valve	Multi-Chamber					
	<p>Arm Splints (Colour blue – blue)</p> <table border="0"> <tr> <td><b>VS62X20V</b></td> <td>Hook &amp; Loop straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS63X20V</b></td> <td>Hook &amp; Loop straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5    Width top/foot: 43-70cm    Length: 56cm    Weight: 0,9kg</p>	<b>VS62X20V</b>	Hook & Loop straps, F-Valve	1-Chamber	<b>VS63X20V</b>	Hook & Loop straps, F-Valve	Multi-Chamber
<b>VS62X20V</b>	Hook & Loop straps, F-Valve	1-Chamber					
<b>VS63X20V</b>	Hook & Loop straps, F-Valve	Multi-Chamber					
	<p>Leg Splints (Colour blue – blue)</p> <table border="0"> <tr> <td><b>VS64X20V</b></td> <td>Hook &amp; Loop straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS65X20V</b></td> <td>Hook &amp; Loop straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5    Width top/foot: 64-70cm    Length: 103cm    Weight: 1,9kg</p>	<b>VS64X20V</b>	Hook & Loop straps, F-Valve	1-Chamber	<b>VS65X20V</b>	Hook & Loop straps, F-Valve	Multi-Chamber
<b>VS64X20V</b>	Hook & Loop straps, F-Valve	1-Chamber					
<b>VS65X20V</b>	Hook & Loop straps, F-Valve	Multi-Chamber					

## MOZART splints with CLICK & replaceable PP or PVC straps

    	<p><b>Forearm splints (X03 = red-blue, X05 = blue-blue)</b></p> <table border="0"> <tr> <td><b>VS60CX03 / VS60CX05</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS60CX03V / VS60CX05V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS60CHX03 / VS60CHX05</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS60CHX05V / VS60CHX05V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td colspan="3"> </td> </tr> <tr> <td><b>VS61CX03 / VS61CX05</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS61CX03V / VS61CX05V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS61CHX03 / VS61CHX05</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS61CHX05V / VS61CHX05V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 3 Width: 35cm Length: 50cm Weight: 0,33kg</p> <p><b>Arm splints (X03 = red-blue, X05 = blue-blue)</b></p> <table border="0"> <tr> <td><b>VS62CX03 / VS62CX05</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS62CX03V / VS62CX05V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS62CHX03 / VS62CHX05</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS62CHX05V / VS62CHX05V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td colspan="3"> </td> </tr> <tr> <td><b>VS63CX03 / VS63CX05</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS63CX03V / VS63CX05V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS63CHX03 / VS63CHX05</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS63CHX05V / VS63CHX05V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width hand/elbow: 43-70cm Length: 56cm Weight: 0,9kg</p> <p><b>Leg Splints (X03 = red-blue, X05 = blue-blue)</b></p> <table border="0"> <tr> <td><b>VS64CX03 / VS64CX05</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS64CX03V / VS64CX05V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS64CHX03 / VS64CHX05</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VS64CHX05V / VS64CHX05V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td colspan="3"> </td> </tr> <tr> <td><b>VS65CX03 / VS65CX05</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS65CX03V / VS65CX05V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS65CHX03 / VS65CHX05</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VS65CHX05V / VS65CHX05V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width top/foot: 64-70cm Length: 103cm Weight: 1,9kg</p>	<b>VS60CX03 / VS60CX05</b>	CLICK, PP-Straps, R-Valve	1-Chamber	<b>VS60CX03V / VS60CX05V</b>	CLICK, PP-Straps, F-Valve	1-Chamber	<b>VS60CHX03 / VS60CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber	<b>VS60CHX05V / VS60CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber	 			<b>VS61CX03 / VS61CX05</b>	CLICK, PP-Straps, R-Valve	Multi-Chamber	<b>VS61CX03V / VS61CX05V</b>	CLICK, PP-Straps, F-Valve	Multi-Chamber	<b>VS61CHX03 / VS61CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VS61CHX05V / VS61CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber	<b>VS62CX03 / VS62CX05</b>	CLICK, PP-Straps, R-Valve	1-Chamber	<b>VS62CX03V / VS62CX05V</b>	CLICK, PP-Straps, F-Valve	1-Chamber	<b>VS62CHX03 / VS62CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber	<b>VS62CHX05V / VS62CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber	 			<b>VS63CX03 / VS63CX05</b>	CLICK, PP-Straps, R-Valve	Multi-Chamber	<b>VS63CX03V / VS63CX05V</b>	CLICK, PP-Straps, F-Valve	Multi-Chamber	<b>VS63CHX03 / VS63CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VS63CHX05V / VS63CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber	<b>VS64CX03 / VS64CX05</b>	CLICK, PP-Straps, R-Valve	1-Chamber	<b>VS64CX03V / VS64CX05V</b>	CLICK, PP-Straps, F-Valve	1-Chamber	<b>VS64CHX03 / VS64CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber	<b>VS64CHX05V / VS64CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber	 			<b>VS65CX03 / VS65CX05</b>	CLICK, PP-Straps, R-Valve	Multi-Chamber	<b>VS65CX03V / VS65CX05V</b>	CLICK, PP-Straps, F-Valve	Multi-Chamber	<b>VS65CHX03 / VS65CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VS65CHX05V / VS65CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber
<b>VS60CX03 / VS60CX05</b>	CLICK, PP-Straps, R-Valve	1-Chamber																																																																																
<b>VS60CX03V / VS60CX05V</b>	CLICK, PP-Straps, F-Valve	1-Chamber																																																																																
<b>VS60CHX03 / VS60CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber																																																																																
<b>VS60CHX05V / VS60CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber																																																																																
<b>VS61CX03 / VS61CX05</b>	CLICK, PP-Straps, R-Valve	Multi-Chamber																																																																																
<b>VS61CX03V / VS61CX05V</b>	CLICK, PP-Straps, F-Valve	Multi-Chamber																																																																																
<b>VS61CHX03 / VS61CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber																																																																																
<b>VS61CHX05V / VS61CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber																																																																																
<b>VS62CX03 / VS62CX05</b>	CLICK, PP-Straps, R-Valve	1-Chamber																																																																																
<b>VS62CX03V / VS62CX05V</b>	CLICK, PP-Straps, F-Valve	1-Chamber																																																																																
<b>VS62CHX03 / VS62CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber																																																																																
<b>VS62CHX05V / VS62CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber																																																																																
<b>VS63CX03 / VS63CX05</b>	CLICK, PP-Straps, R-Valve	Multi-Chamber																																																																																
<b>VS63CX03V / VS63CX05V</b>	CLICK, PP-Straps, F-Valve	Multi-Chamber																																																																																
<b>VS63CHX03 / VS63CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber																																																																																
<b>VS63CHX05V / VS63CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber																																																																																
<b>VS64CX03 / VS64CX05</b>	CLICK, PP-Straps, R-Valve	1-Chamber																																																																																
<b>VS64CX03V / VS64CX05V</b>	CLICK, PP-Straps, F-Valve	1-Chamber																																																																																
<b>VS64CHX03 / VS64CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	1-Chamber																																																																																
<b>VS64CHX05V / VS64CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	1-Chamber																																																																																
<b>VS65CX03 / VS65CX05</b>	CLICK, PP-Straps, R-Valve	Multi-Chamber																																																																																
<b>VS65CX03V / VS65CX05V</b>	CLICK, PP-Straps, F-Valve	Multi-Chamber																																																																																
<b>VS65CHX03 / VS65CHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber																																																																																
<b>VS65CHX05V / VS65CHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber																																																																																
    	<p><b>Forearm splints (X03 = red-blue, X05 = blue-blue)</b></p> <table border="0"> <tr> <td><b>VS61CPX03 / VS61CPX05</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VS61CPX03V / VS61CPX05V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VS61CPHX03 / VS61CPHX05</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VS61CPHX05V / VS61CPHX05V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> </table> <p>Restraint straps: 3 Width: 35cm Length: 50cm Weight: 0,33kg</p> <p><b>Arm splints (X03 = red-blue, X05 = blue-blue)</b></p> <table border="0"> <tr> <td><b>VS63CPX03 / VS63CPX05</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VS63CPX03V / VS63CPX05V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VS63CPHX03 / VS63CPHX05</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VS63CPHX05V / VS63CPHX05V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width hand/elbow: 43-70cm Length: 56cm Weight: 0,9kg</p> <p><b>Leg Splints (X03 = red-blue, X05 = blue-blue)</b></p> <table border="0"> <tr> <td><b>VS65CPX03 / VS65CPX05</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VS65CPX03V / VS65CPX05V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VS65CPHX03 / VS65CPHX05</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VS65CPHX05V / VS65CPHX05V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width top/foot: 64-70cm Length: 103cm Weight: 1,9kg</p>	<b>VS61CPX03 / VS61CPX05</b>	CLICK, PP-Straps, R-Valve	DOT-Chamber	<b>VS61CPX03V / VS61CPX05V</b>	CLICK, PP-Straps, F-Valve	DOT-Chamber	<b>VS61CPHX03 / VS61CPHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber	<b>VS61CPHX05V / VS61CPHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber	<b>VS63CPX03 / VS63CPX05</b>	CLICK, PP-Straps, R-Valve	DOT-Chamber	<b>VS63CPX03V / VS63CPX05V</b>	CLICK, PP-Straps, F-Valve	DOT-Chamber	<b>VS63CPHX03 / VS63CPHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber	<b>VS63CPHX05V / VS63CPHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber	<b>VS65CPX03 / VS65CPX05</b>	CLICK, PP-Straps, R-Valve	DOT-Chamber	<b>VS65CPX03V / VS65CPX05V</b>	CLICK, PP-Straps, F-Valve	DOT-Chamber	<b>VS65CPHX03 / VS65CPHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber	<b>VS65CPHX05V / VS65CPHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																																													
<b>VS61CPX03 / VS61CPX05</b>	CLICK, PP-Straps, R-Valve	DOT-Chamber																																																																																
<b>VS61CPX03V / VS61CPX05V</b>	CLICK, PP-Straps, F-Valve	DOT-Chamber																																																																																
<b>VS61CPHX03 / VS61CPHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber																																																																																
<b>VS61CPHX05V / VS61CPHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																																																																																
<b>VS63CPX03 / VS63CPX05</b>	CLICK, PP-Straps, R-Valve	DOT-Chamber																																																																																
<b>VS63CPX03V / VS63CPX05V</b>	CLICK, PP-Straps, F-Valve	DOT-Chamber																																																																																
<b>VS63CPHX03 / VS63CPHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber																																																																																
<b>VS63CPHX05V / VS63CPHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																																																																																
<b>VS65CPX03 / VS65CPX05</b>	CLICK, PP-Straps, R-Valve	DOT-Chamber																																																																																
<b>VS65CPX03V / VS65CPX05V</b>	CLICK, PP-Straps, F-Valve	DOT-Chamber																																																																																
<b>VS65CPHX03 / VS65CPHX05</b>	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber																																																																																
<b>VS65CPHX05V / VS65CPHX05V</b>	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																																																																																





	<p><b>Arm Splints (Colour blue – blue)</b></p> <table border="0"> <tbody> <tr> <td><b>VA14X02</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VA14X02V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VA14HX02</b></td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VA14HX02V</b></td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VA14X01</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14X01V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14HX01</b></td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14HX01V</b></td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </tbody> </table> <p>Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg</p>	<b>VA14X02</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber	<b>VA14X02V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber	<b>VA14HX02</b>	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber	<b>VA14HX02V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber	<b>VA14X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	<b>VA14X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber	<b>VA14HX01</b>	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VA14HX01V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber
<b>VA14X02</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber																							
<b>VA14X02V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber																							
<b>VA14HX02</b>	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber																							
<b>VA14HX02V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber																							
<b>VA14X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber																							
<b>VA14X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber																							
<b>VA14HX01</b>	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
<b>VA14HX01V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber																							
	<p><b>Ankle Splints (Colour blue – blue)</b></p> <table border="0"> <tbody> <tr> <td><b>VK15X02</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VK15X02V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VK15HX02</b></td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VK15HX02V</b></td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VK15X01</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VK15X01V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA15HX01</b></td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA15HX01V</b></td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </tbody> </table> <p>Restraint straps: 4 Width: 47cm Length: 57cm Weight: 0,9kg</p>	<b>VK15X02</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber	<b>VK15X02V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber	<b>VK15HX02</b>	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber	<b>VK15HX02V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber	<b>VK15X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	<b>VK15X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber	<b>VA15HX01</b>	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VA15HX01V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber
<b>VK15X02</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber																							
<b>VK15X02V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber																							
<b>VK15HX02</b>	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber																							
<b>VK15HX02V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber																							
<b>VK15X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber																							
<b>VK15X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber																							
<b>VA15HX01</b>	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
<b>VA15HX01V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber																							
	<p><b>Leg Splints</b></p> <table border="0"> <tbody> <tr> <td><b>VL21X02</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL21X02V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL21HX02</b></td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL21HX02V</b></td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL21X01</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21X01V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21HX01</b></td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21HX01V</b></td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </tbody> </table> <p>Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg</p>	<b>VL21X02</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber	<b>VL21X02V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber	<b>VL21HX02</b>	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber	<b>VL21HX02V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber	<b>VL21X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	<b>VL21X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber	<b>VL21HX01</b>	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VL21HX01V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber
<b>VL21X02</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber																							
<b>VL21X02V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber																							
<b>VL21HX02</b>	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber																							
<b>VL21HX02V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber																							
<b>VL21X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber																							
<b>VL21X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber																							
<b>VL21HX01</b>	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
<b>VL21HX01V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber																							
	<p><b>Large Leg Splints (Colour blue – blue)</b></p> <table border="0"> <tbody> <tr> <td><b>VL18X02</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL18X02V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL18HX02</b></td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL18HX02V</b></td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL18X01</b></td> <td>Red Stripes, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL18X01V</b></td> <td>Red Stripes, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL18HX01</b></td> <td>Red Stripes, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL18HX01V</b></td> <td>Red Stripes, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </tbody> </table> <p>Restraint straps: 5 Width: top/hand: 70-45cm Length: 135cm Weight: 2,1kg</p>	<b>VL18X02</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber	<b>VL18X02V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber	<b>VL18HX02</b>	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber	<b>VL18HX02V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber	<b>VL18X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber	<b>VL18X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber	<b>VL18HX01</b>	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VL18HX01V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber
<b>VL18X02</b>	Red Stripes, PP-Straps, R-Valve	1-Chamber																							
<b>VL18X02V</b>	Red Stripes, PP-Straps, F-Valve	1-Chamber																							
<b>VL18HX02</b>	Red Stripes, Wipe-clean-Straps, R-Valve	1-Chamber																							
<b>VL18HX02V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	1-Chamber																							
<b>VL18X01</b>	Red Stripes, PP-Straps, R-Valve	Multi-Chamber																							
<b>VL18X01V</b>	Red Stripes, PP-Straps, F-Valve	Multi-Chamber																							
<b>VL18HX01</b>	Red Stripes, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
<b>VL18HX01V</b>	Red Stripes, Wipe-clean-Straps, F-Valve	Multi-Chamber																							



## VIKING splints with Cleats & replaceable PP or PVC straps

          	<p><b>Arm Splints (Colour blue – blue)</b></p> <table border="0"> <tr> <td><b>VA14X04</b></td> <td>CLEATS, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VA14X04V</b></td> <td>CLEATS, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VA14HX04</b></td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VA14HX04V</b></td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VA14X03</b></td> <td>CLEATS, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14X03V</b></td> <td>CLEATS, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14HX03</b></td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14HX03V</b></td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg</p> <p><b>Leg Splints (Colour blue – blue)</b></p> <table border="0"> <tr> <td><b>VL21X04</b></td> <td>CLEATS, PP-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL21X04V</b></td> <td>CLEATS, PP-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL21HX04</b></td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL21HX04V</b></td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>1-Chamber</td> </tr> <tr> <td><b>VL21X03</b></td> <td>CLEATS, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21X03V</b></td> <td>CLEATS, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21HX31</b></td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21HX03V</b></td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg</p>	<b>VA14X04</b>	CLEATS, PP-Straps, R-Valve	1-Chamber	<b>VA14X04V</b>	CLEATS, PP-Straps, F-Valve	1-Chamber	<b>VA14HX04</b>	CLEATS, Wipe-clean-Straps, R-Valve	1-Chamber	<b>VA14HX04V</b>	CLEATS, Wipe-clean-Straps, F-Valve	1-Chamber	<b>VA14X03</b>	CLEATS, PP-Straps, R-Valve	Multi-Chamber	<b>VA14X03V</b>	CLEATS, PP-Straps, F-Valve	Multi-Chamber	<b>VA14HX03</b>	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VA14HX03V</b>	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber	<b>VL21X04</b>	CLEATS, PP-Straps, R-Valve	1-Chamber	<b>VL21X04V</b>	CLEATS, PP-Straps, F-Valve	1-Chamber	<b>VL21HX04</b>	CLEATS, Wipe-clean-Straps, R-Valve	1-Chamber	<b>VL21HX04V</b>	CLEATS, Wipe-clean-Straps, F-Valve	1-Chamber	<b>VL21X03</b>	CLEATS, PP-Straps, R-Valve	Multi-Chamber	<b>VL21X03V</b>	CLEATS, PP-Straps, F-Valve	Multi-Chamber	<b>VL21HX31</b>	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VL21HX03V</b>	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber
<b>VA14X04</b>	CLEATS, PP-Straps, R-Valve	1-Chamber																																															
<b>VA14X04V</b>	CLEATS, PP-Straps, F-Valve	1-Chamber																																															
<b>VA14HX04</b>	CLEATS, Wipe-clean-Straps, R-Valve	1-Chamber																																															
<b>VA14HX04V</b>	CLEATS, Wipe-clean-Straps, F-Valve	1-Chamber																																															
<b>VA14X03</b>	CLEATS, PP-Straps, R-Valve	Multi-Chamber																																															
<b>VA14X03V</b>	CLEATS, PP-Straps, F-Valve	Multi-Chamber																																															
<b>VA14HX03</b>	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber																																															
<b>VA14HX03V</b>	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber																																															
<b>VL21X04</b>	CLEATS, PP-Straps, R-Valve	1-Chamber																																															
<b>VL21X04V</b>	CLEATS, PP-Straps, F-Valve	1-Chamber																																															
<b>VL21HX04</b>	CLEATS, Wipe-clean-Straps, R-Valve	1-Chamber																																															
<b>VL21HX04V</b>	CLEATS, Wipe-clean-Straps, F-Valve	1-Chamber																																															
<b>VL21X03</b>	CLEATS, PP-Straps, R-Valve	Multi-Chamber																																															
<b>VL21X03V</b>	CLEATS, PP-Straps, F-Valve	Multi-Chamber																																															
<b>VL21HX31</b>	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber																																															
<b>VL21HX03V</b>	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber																																															
          	<p><b>Arm Splints (Colour red – blue)</b></p> <table border="0"> <tr> <td><b>VA14X07</b></td> <td>CLEATS, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14X07V</b></td> <td>CLEATS, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14HX07</b></td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14HX07V</b></td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg</p> <p><b>Leg Splints (Colour red – blue)</b></p> <table border="0"> <tr> <td><b>VL21X07</b></td> <td>CLEATS, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21X07V</b></td> <td>CLEATS, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21HX07</b></td> <td>CLEATS, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21HX03V</b></td> <td>CLEATS, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg</p>	<b>VA14X07</b>	CLEATS, PP-Straps, R-Valve	Multi-Chamber	<b>VA14X07V</b>	CLEATS, PP-Straps, F-Valve	Multi-Chamber	<b>VA14HX07</b>	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VA14HX07V</b>	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber	<b>VL21X07</b>	CLEATS, PP-Straps, R-Valve	Multi-Chamber	<b>VL21X07V</b>	CLEATS, PP-Straps, F-Valve	Multi-Chamber	<b>VL21HX07</b>	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VL21HX03V</b>	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber																								
<b>VA14X07</b>	CLEATS, PP-Straps, R-Valve	Multi-Chamber																																															
<b>VA14X07V</b>	CLEATS, PP-Straps, F-Valve	Multi-Chamber																																															
<b>VA14HX07</b>	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber																																															
<b>VA14HX07V</b>	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber																																															
<b>VL21X07</b>	CLEATS, PP-Straps, R-Valve	Multi-Chamber																																															
<b>VL21X07V</b>	CLEATS, PP-Straps, F-Valve	Multi-Chamber																																															
<b>VL21HX07</b>	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber																																															
<b>VL21HX03V</b>	CLEATS, Wipe-clean-Straps, F-Valve	Multi-Chamber																																															




## VIKING splints with CLICK & replaceable PP or PVC straps

          	<p><b>Arm Splints (X07 = red-blue, X03 = blue-blue)</b></p> <table border="0"> <tr> <td><b>VA14CX07 / VA14CX03</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14CX07V / VA14CX03V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14CHX07 / VA14CHX03</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VA14CHX07V / VA14CHX03V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg</p> <p><b>Leg Splints (X07 = red-blue, X03 = blue-blue)</b></p> <table border="0"> <tr> <td><b>VL21CX07 / VL21CX03</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21CX07V / VL21CX03V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21CHX07 / VL21CHX03</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>Multi-Chamber</td> </tr> <tr> <td><b>VL21CHX07V / VL21CHX03V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>Multi-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg</p>	<b>VA14CX07 / VA14CX03</b>	CLICK, PP-Straps, R-Valve	Multi-Chamber	<b>VA14CX07V / VA14CX03V</b>	CLICK, PP-Straps, F-Valve	Multi-Chamber	<b>VA14CHX07 / VA14CHX03</b>	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VA14CHX07V / VA14CHX03V</b>	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber	<b>VL21CX07 / VL21CX03</b>	CLICK, PP-Straps, R-Valve	Multi-Chamber	<b>VL21CX07V / VL21CX03V</b>	CLICK, PP-Straps, F-Valve	Multi-Chamber	<b>VL21CHX07 / VL21CHX03</b>	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber	<b>VL21CHX07V / VL21CHX03V</b>	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber
<b>VA14CX07 / VA14CX03</b>	CLICK, PP-Straps, R-Valve	Multi-Chamber																							
<b>VA14CX07V / VA14CX03V</b>	CLICK, PP-Straps, F-Valve	Multi-Chamber																							
<b>VA14CHX07 / VA14CHX03</b>	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
<b>VA14CHX07V / VA14CHX03V</b>	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber																							
<b>VL21CX07 / VL21CX03</b>	CLICK, PP-Straps, R-Valve	Multi-Chamber																							
<b>VL21CX07V / VL21CX03V</b>	CLICK, PP-Straps, F-Valve	Multi-Chamber																							
<b>VL21CHX07 / VL21CHX03</b>	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber																							
<b>VL21CHX07V / VL21CHX03V</b>	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber																							
          	<p><b>Arm Splints (X07 = red-blue, X03 = blue-blue)</b></p> <table border="0"> <tr> <td><b>VA14CPX07 / VA14CPX03</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VA14CPX07V / VA14CPX03V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VA14CPHX07 / VA14CPHX03</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VA14CPHX07V / VA14CPHX03V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> </table> <p>Restraint straps: 4 Width: top/hand: 50-30cm Length: 67cm Weight: 0,9kg</p> <p><b>Leg Splints</b></p> <table border="0"> <tr> <td><b>VL21CPX07 / VL21CPX03</b></td> <td>CLICK, PP-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VL21CPX07V / VL21CPX03V</b></td> <td>CLICK, PP-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VL21CPHX07 / VL21CPHX03</b></td> <td>CLICK, Wipe-clean-Straps, R-Valve</td> <td>DOT-Chamber</td> </tr> <tr> <td><b>VL21CPHX07V / VL21CPHX03V</b></td> <td>CLICK, Wipe-clean-Straps, F-Valve</td> <td>DOT-Chamber</td> </tr> </table> <p>Restraint straps: 5 Width: top/hand: 70-40cm Length: 110cm Weight: 1,9kg</p>	<b>VA14CPX07 / VA14CPX03</b>	CLICK, PP-Straps, R-Valve	DOT-Chamber	<b>VA14CPX07V / VA14CPX03V</b>	CLICK, PP-Straps, F-Valve	DOT-Chamber	<b>VA14CPHX07 / VA14CPHX03</b>	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber	<b>VA14CPHX07V / VA14CPHX03V</b>	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber	<b>VL21CPX07 / VL21CPX03</b>	CLICK, PP-Straps, R-Valve	DOT-Chamber	<b>VL21CPX07V / VL21CPX03V</b>	CLICK, PP-Straps, F-Valve	DOT-Chamber	<b>VL21CPHX07 / VL21CPHX03</b>	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber	<b>VL21CPHX07V / VL21CPHX03V</b>	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber
<b>VA14CPX07 / VA14CPX03</b>	CLICK, PP-Straps, R-Valve	DOT-Chamber																							
<b>VA14CPX07V / VA14CPX03V</b>	CLICK, PP-Straps, F-Valve	DOT-Chamber																							
<b>VA14CPHX07 / VA14CPHX03</b>	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber																							
<b>VA14CPHX07V / VA14CPHX03V</b>	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																							
<b>VL21CPX07 / VL21CPX03</b>	CLICK, PP-Straps, R-Valve	DOT-Chamber																							
<b>VL21CPX07V / VL21CPX03V</b>	CLICK, PP-Straps, F-Valve	DOT-Chamber																							
<b>VL21CPHX07 / VL21CPHX03</b>	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber																							
<b>VL21CPHX07V / VL21CPHX03V</b>	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber																							




## RedVac Vacuum Splints "GLOBAL"

### GLOBAL splints with CLEATS & replaceable Loop straps




7

    	<b>Forearm Splints (Colour blue – blue)</b>		
	<b>VS100X02</b>	CLEATS, Loop straps, R-Valve	1-Chamber
	<b>VS100X02V</b>	CLEATS, Loop straps, F-Valve	1-Chamber
	<b>VS100X01</b>	CLEATS, Loop straps, R-Valve	Multi-Chamber
	<b>VS100X01V</b>	CLEATS, Loop straps, F-Valve	Multi-Chamber
	Restraint straps: 2 Width: 35cm Length: 53cm Weight: 0,5kg		
	<b>Arm Splints (Colour blue – blue)</b>		
	<b>VS120X02</b>	CLEATS, Loop straps, R-Valve	1-Chamber
	<b>VS120X02V</b>	CLEATS, Loop straps, F-Valve	1-Chamber
	<b>VS120X01</b>	CLEATS, Loop straps, R-Valve	Multi-Chamber
	<b>VS120X01V</b>	CLEATS, Loop straps, F-Valve	Multi-Chamber
	Restraint straps: 4 Width: 63-47cm Length: 74cm Weight: 1,1kg		
<b>Leg Splints (Colour blue – blue)</b>			
<b>VS140X02</b>	CLEATS, Loop straps, R-Valve	1-Chamber	
<b>VS140X02V</b>	CLEATS, Loop straps, F-Valve	1-Chamber	
<b>VS140X01</b>	CLEATS, Loop straps, R-Valve	Multi-Chamber	
<b>VS140X01V</b>	CLEATS, Loop straps, F-Valve	Multi-Chamber	
Restraint straps: 6 Width: 70-36cm Length: 110cm Weight: 2,1kg			




### GLOBAL splints with Red Stripes & replaceable Loop straps

    	<b>Forearm Splints</b>		
	<b>VS100X04</b>	Red Stripes, Loop straps, R-Valve	1-Chamber
	<b>VS100X04V</b>	Red Stripes, Loop straps, F-Valve	1-Chamber
	<b>VS100X03</b>	Red Stripes, Loop straps, R-Valve	Multi-Chamber
	<b>VS100X03V</b>	Red Stripes, Loop straps, R-Valve	Multi-Chamber
	Restraint straps: 2 Width: 35cm Length: 53cm Weight: 0,5kg		
	<b>Arm Splints</b>		
	<b>VS120X04</b>	Red Stripes, Loop straps, R-Valve	1-Chamber
	<b>VS120X04V</b>	Red Stripes, Loop straps, F-Valve	1-Chamber
	<b>VS120X03</b>	Red Stripes, Loop straps, R-Valve	Multi-Chamber
	<b>VS120X03V</b>	Red Stripes, Loop straps, R-Valve	Multi-Chamber
	Restraint straps: 4 Width: 63-47cm Length: 74cm Weight: 1,1kg		
<b>Leg Splints</b>			
<b>VS140X04</b>	Red Stripes, Loop straps, R-Valve	1-Chamber	
<b>VS140X04V</b>	Red Stripes, Loop straps, F-Valve	1-Chamber	
<b>VS140X03</b>	Red Stripes, Loop straps, R-Valve	Multi-Chamber	
<b>VS140X03V</b>	Red Stripes, Loop straps, R-Valve	Multi-Chamber	
Restraint straps: 6 Width: 70-36cm Length: 110cm Weight: 2,1kg			

## GLOBAL splints with CLICK & replaceable Loop straps straps

<p><b>Forearm Splints</b></p> 	<p> <b>VS100CKX04</b> CLICK, Loop straps, R-Valve 1-Chamber  <b>VS100CKX04V</b> CLICK, Loop straps, F-Valve 1-Chamber  <b>VS100CKX03</b> CLICK, Loop straps, R-Valve Multi-Chamber  <b>VS100CKX03V</b> CLICK, Loop straps, F-Valve Multi-Chamber  <b>VS100CPKX03</b> CLICK, Loop straps, R-Valve DOT-Chamber  <b>VS100CPKX03V</b> CLICK, Loop straps, R-Valve DOT-Chamber         </p> <p>Restraint straps: 2 Width: 35cm Length: 53cm Weight: 0,5kg</p>
<p><b>Arm Splints</b></p> 	<p> <b>VS120CKX04</b> CLICK, Loop straps, R-Valve 1-Chamber  <b>VS120CKX04V</b> CLICK, Loop straps, F-Valve 1-Chamber  <b>VS120CKX03</b> CLICK, Loop straps, R-Valve Multi-Chamber  <b>VS120CKX03V</b> CLICK, Loop straps, F-Valve Multi-Chamber  <b>VS120CPKX03</b> CLICK, Loop straps, R-Valve DOT-Chamber  <b>VS120CPKX03V</b> CLICK, Loop straps, R-Valve DOT-Chamber         </p> <p>Restraint straps: 4 Width: 63-47cm Length: 74cm Weight: 1,1kg</p>
<p><b>Leg Splints</b></p> 	<p> <b>VS140CKX04</b> CLICK, Loop straps, R-Valve 1-Chamber  <b>VS140CKX04V</b> CLICK, Loop straps, F-Valve 1-Chamber  <b>VS140CKX03</b> CLICK, Loop straps, R-Valve Multi-Chamber  <b>VS140CKX03V</b> CLICK, Loop straps, F-Valve Multi-Chamber  <b>VS140CPKX03</b> CLICK, Loop straps, R-Valve DOT-Chamber  <b>VS140CPKX03V</b> CLICK, Loop straps, R-Valve DOT-Chamber         </p> <p>Restraint straps: 6 Width: 70-36cm Length: 110cm Weight: 2,1kg</p>


## GLOBAL splints with CLICK & replaceable PP or PVC straps

<p><b>Forearm Splints</b></p> 	<p> <b>VS100CX04</b> CLICK, PP-Straps, R-Valve 1-Chamber  <b>VS100CX04V</b> CLICK, PP-Straps, F-Valve 1-Chamber  <b>VS100CHX04</b> CLICK, Wipe-clean-Straps, R-Valve 1-Chamber  <b>VS100CHX04V</b> CLICK, Wipe-clean-Straps, F-Valve 1-Chamber  <b>VS100CX03</b> CLICK, PP-Straps, R-Valve Multi-Chamber  <b>VS100CX03V</b> CLICK, PP-Straps, F-Valve Multi-Chamber  <b>VS100CHX03</b> CLICK, Wipe-clean-Straps, R-Valve Multi-Chamber  <b>VS100CHX03V</b> CLICK, Wipe-clean-Straps, F-Valve Multi-Chamber  <b>VS100CPX03</b> CLICK, PP-Straps, R-Valve DOT-Chamber  <b>VS100CPX03V</b> CLICK, PP-Straps, F-Valve DOT-Chamber  <b>VS100CPHX03</b> CLICK, Wipe-clean-Straps, R-Valve DOT-Chamber  <b>VS100CPHX03V</b> CLICK, Wipe-clean-Straps, F-Valve DOT-Chamber         </p> <p>Restraint straps: 3 Width: 35cm Length: 53cm Weight: 0,5kg</p>
<p><b>Arm Splints</b></p> 	<p> <b>VS120CX04</b> CLICK, PP-Straps, R-Valve 1-Chamber  <b>VS120CX04V</b> CLICK, PP-Straps, F-Valve 1-Chamber  <b>VS120CHX04</b> CLICK, Wipe-clean-Straps, R-Valve 1-Chamber  <b>VS120CHX04V</b> CLICK, Wipe-clean-Straps, F-Valve 1-Chamber  <b>VS120CX03</b> CLICK, PP-Straps, R-Valve Multi-Chamber  <b>VS120CX03V</b> CLICK, PP-Straps, F-Valve Multi-Chamber  <b>VS120CHX03</b> CLICK, Wipe-clean-Straps, R-Valve Multi-Chamber  <b>VS120CHX03V</b> CLICK, Wipe-clean-Straps, F-Valve Multi-Chamber  <b>VS120CPX03</b> CLICK, PP-Straps, R-Valve DOT-Chamber  <b>VS120CPX03V</b> CLICK, PP-Straps, F-Valve DOT-Chamber  <b>VS120CPHX03</b> CLICK, Wipe-clean-Straps, R-Valve DOT-Chamber  <b>VS120CPHX03V</b> CLICK, Wipe-clean-Straps, F-Valve DOT-Chamber         </p> <p>Restraint straps: 4 Width: 63-47cm Length: 74cm Weight: 1,1kg</p>
<p><b>Leg Splints</b></p> 	<p> <b>VS140CX04</b> CLICK, PP-Straps, R-Valve 1-Chamber  <b>VS140CX04V</b> CLICK, PP-Straps, F-Valve 1-Chamber  <b>VS140CHX04</b> CLICK, Wipe-clean-Straps, R-Valve 1-Chamber  <b>VS140CHX04V</b> CLICK, Wipe-clean-Straps, F-Valve 1-Chamber  <b>VS140CX03</b> CLICK, PP-Straps, R-Valve Multi-Chamber  <b>VS140CX03V</b> CLICK, PP-Straps, F-Valve Multi-Chamber  <b>VS140CHX03</b> CLICK, Wipe-clean-Straps, R-Valve Multi-Chamber  <b>VS140CHX03V</b> CLICK, Wipe-clean-Straps, F-Valve Multi-Chamber  <b>VS140CPX03</b> CLICK, PP-Straps, R-Valve DOT-Chamber  <b>VS140CPX03V</b> CLICK, PP-Straps, F-Valve DOT-Chamber  <b>VS140CPHX03</b> CLICK, Wipe-clean-Straps, R-Valve DOT-Chamber  <b>VS140CPHX03V</b> CLICK, Wipe-clean-Straps, F-Valve DOT-Chamber         </p> <p>Restraint straps: 6 Width: 70-36cm Length: 110cm Weight: 2,1kg</p>



## RedVac Vacuum Splints “REMBRANDT”

### REMBRANDT splints with Loop straps



	<p><b>Arm Splint</b></p> <p><b>VS101X02</b> Hook &amp; Loop straps, R-Valve Multi-Chamber</p> <p>Restraint straps: 2            Width: Top/Hand: 40-30cm            Length: 73cm            Weight: 0,7kg</p>
	<p><b>Leg Splint</b></p> <p><b>VS103X02</b> Hook &amp; Loop straps, R-Valve Multi-Chamber</p> <p>Restraint straps: 4            Width: Top/Foot: 60-18cm            Length: 96cm            Weight: 1,9kg</p>
	<p><b>Leg Splint</b></p> <p><b>VS105X02</b> Hook &amp; Loop straps, R-Valve Multi-Chamber</p> <p>Restraint straps: 5            Width: Top/Foot: 73-18cm            Length: 133cm            Weight: 2,2kg</p>

## RedVac Vacuum Splints “BASIC”



	<p><b>Arm Splint</b></p> <p><b>VSDS</b> Basic, without straps, F-Valve            Width: Top/Hand: 40-30cm Length: 73cm Weight: 0,7kg</p> <p><b>VSDSX01</b> Basic with Hook &amp; Loop straps, F-Valve            Restraint straps: 3            Width: Top/Hand: 40-30cm Length: 73cm Weight: 0,7kg</p>
	<p><b>Leg Splint</b></p> <p><b>VSDM</b> Basic, without straps, F-Valve            Width: Top/Foot: 60-18cm Length: 96cm Weight: 1,9kg</p> <p><b>VSDMX01</b> Basic with Hook &amp; Loop straps, F-Valve            Restraint straps: 4            Width: Top/Foot: 60-18cm Length: 96cm Weight: 1,9kg</p>
	<p><b>Leg Splint</b></p> <p><b>VSDL</b> Basic, without straps, F-Valve            Width: Top/Foot: 73-18cm Length: 133cm Weight: 2,2kg</p> <p><b>VSDLX01</b> Basic with Hook &amp; Loop straps, F-Valve            Restraint straps: 5            Width: Top/Foot: 73-18cm Length: 133cm Weight: 2,2kg</p>

# RedVac Pediatric Vacuum Mattresses

## Vacuum Mattresses with CLEATS & PP or PVC straps



	VM710X22	CLEATS, PP-Straps, R-Valve	1-Chamber	
	VM710X22V	CLEATS, PP-Straps, F-Valve	1-Chamber	
	VM710HX22	CLEATS, Wipe-clean-Straps, R-Valve	1-Chamber	
	VM710XH22V	CLEATS, Wipe-clean-Straps, F-Valv	1-Chamber	
	VM710X08	CLEATS, PP-Straps, R-Valve	Multi-Chamber	
	VM710X08V	CLEATS, PP-Straps, F-Valve	Multi-Chamber	
	VM710HX08	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber	
	VM710XH08V	CLEATS, Wipe-clean-Straps, F-Valv	Multi-Chamber	
Carry handles: 4 Restraint straps: 3 Width: 70cm Length: 100cm Weight: 3,5kg				
	VM710X23	CLEATS, PP-Straps, R-Valve	Multi-Chamber	
	VM710X23V	CLEATS, PP-Straps, F-Valve	Multi-Chamber	
	VM710HX23	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber	
	VM710HX23V	CLEATS, Wipe-clean-Straps, F-Valv	Multi-Chamber	
	Carry handles: 4 Restraint straps: 3 Width: 70cm Length: 100cm Weight: 3,5kg			
	VM712X09	CLEATS, PP-Straps, R-Valve	Multi-Chamber	
	VM712X09V	CLEATS, PP-Straps, F-Valve	Multi-Chamber	
	VM712HX09	CLEATS, Wipe-clean-Straps, R-Valve	Multi-Chamber	
VM712HX09V	CLEATS, Wipe-clean-Straps, F-Valv	Multi-Chamber		
Carry handles: 4 Restraint straps: 3 Width: 70cm Length: 100cm Weight: 3,5kg				

## Vacuum Mattresses with CLICK & PP or PVC straps

	VM710CX23	CLICK, PP-Straps, R-Valve	Multi-Chamber
	VM710CX23V	CLICK, PP-Straps, F-Valve	Multi-Chamber
	VM710CHX23	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber
	VM710CHX23V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber
	VM710CPX23	CLICK, PP-Straps, R-Valve	DOT-Chamber
	VM710CPX23V	CLICK, PP-Straps, F-Valve	DOT-Chamber
	VM710CPHX23	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber
	VM710CPHX23V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber
Carry handles: 4 Restraint straps: 3 Width: 70cm Length: 100cm Weight: 3,5kg			
	VM712CX23	CLICK, PP-Straps, R-Valve	Multi-Chamber
	VM712CX23V	CLICK, PP-Straps, F-Valve	Multi-Chamber
	VM712CHX23	CLICK, Wipe-clean-Straps, R-Valve	Multi-Chamber
	VM712CHX23V	CLICK, Wipe-clean-Straps, F-Valve	Multi-Chamber
	VM712CPX23	CLICK, PP-Straps, R-Valve	DOT-Chamber
	VM712CPX23V	CLICK, PP-Straps, F-Valve	DOT-Chamber
	VM712CPHX23	CLICK, Wipe-clean-Straps, R-Valve	DOT-Chamber
	VM712CPHX23V	CLICK, Wipe-clean-Straps, F-Valve	DOT-Chamber
Carry handles: 4 Restraint straps: 3 Width: 70cm Length: 120cm Weight: 3,5kg			

# RedVac Pediatric Vacuum Splints

## Splints with replaceable Loop straps

Arm splint, angled 	VS56X20V	Hook & Loop straps, F-Valve	1-Chamber
Restraint straps: 3 Width top/hand: 50-28cm Length: 42cm Weight: 0,6 kg			
Leg Splint 	VS58X20V	Hook & Loop straps, F-Valve	1-Chamber
Restraint straps: 3 Width top/foot: 38-25cm Length: 56cm Weight: 0,7kg			

## Splints with RED STRIPES & replaceable PP or PVC straps

Arm Splint, angled 	VS56X01	Red Stripes, PP-Straps, R-Valve	1-Chamber
	VS56X01V	Red Stripes, PP-Straps, F-Valve	1-Chamber
Restraint straps: 3 Width top/hand: 50-28cm Length: 42cm Weight: 0,6 kg			
Leg Splint 	VS58X02	Red Stripes, PP-Straps, R-Valve	1-Chamber
	VS58X02V	Red Stripes, PP-Straps, F-Valve	1-Chamber
	VS59X01	Red Stripes, PP-Straps, R-Valve	Multi-Chamber
	VS59X01V	Red Stripes, PP-Straps, F-Valve	Multi-Chamber
Restraint straps: 3 Width top/hand: 50-28cm Length: 42cm Weight: 0,6 kg			




# RedVac Pediatric Restraint System for the ambulance gurney

## RedVac JUNIOR

Pediatric Restraint System for use with or without vacuum mattress




**JUNIOR is a restraint system for safe transportation of children in ambulances. It is compatible with most ambulance gurneys. The system is suitable for new-borns and children up to around 7 years with body length 130 cm and weight of 35 kg.**

<p><b>Junior</b></p> 	<p><b>KF01X07</b></p> <p>Width head: 39 cm Width foot: 20 cm Length: 74 cm Weight: 1,3 kg</p> <p>Restraint straps: 4 Securing straps: 6</p>		
<p><b>Vacuum mattresses for JUNIOR</b></p> 	<p><b>VKF01PX02</b> Hook &amp; Loop straps, R-Valve DOT-Chamber  <b>VKF01PX02V</b> Hook &amp; Loop straps, F-Valve DOT-Chamber  <b>VKF01CPX02</b> CLICK with PP-Straps, R-Valve DOT-Chamber  <b>VKF01CPX02V</b> CLICK with PP-Straps, F-Valve DOT-Chamber  <b>VKF01CPHX02</b> CLICK with Wipe-clean-Straps, R-Valve, R-Valve DOT-Chamber  <b>VKF01CPHX02V</b> CLICK with Wipe-clean-Straps, F-Valve DOT-Chamber</p> <p>Width: 45 cm Length: 80 cm Weight: 0,75 kg Restraint straps: 3</p>		
	<p><b>VKF02PX01</b> Hook &amp; Loop straps, R-Valve DOT-Chamber  <b>VKF02PX01V</b> Hook &amp; Loop straps, F-Valve DOT-Chamber  <b>VKF02CPX01</b> CLICK with PP-Straps, R-Valve DOT-Chamber  <b>VKF02CPX01V</b> CLICK with PP-Straps, F-Valve DOT-Chamber  <b>VKF02CPHX01</b> CLICK with Wipe-clean-Straps, R-Valve DOT-Chamber  <b>VKF02CPHX01V</b> CLICK with Wipe-clean-Straps, F-Valve DOT-Chamber</p> <p>Width Head/Middle-Foot: 45/62 cm Length: 138 cm Weight: 1,2 kg Restraint straps: 5</p>		

## Vacuum Pumps & Carry Bags

	<p><b>HP01</b> One-Hand pump, plastic Length: 21cm Width: 14cm Weight: 0,6kg</p>
	<p><b>HP03</b> Hand pump, aluminium Length: 30cm Width: 4cm Weight: 0,4kg</p>
	<p><b>HP04</b> Hand pump, plastic Length: 28cm Width: 11 / 4,5cm Weight: 0,170kg</p>
	<p><b>FP01</b> Foot pump, plastic Length: 43cm Width: 11cm Weight: 0,75kg</p>
	<p><b>HFP02</b> Hand- and Foot pump, plastic. Length: 45cm Width: 6cm Weight: 0,95kg</p>
	<p><b>FP04</b> Foot pump, aluminium Length: 54cm Width: 6cm Weight: 1,4kg</p>
<p>Carry bag for RedVac Junior</p> 	<p><b>TS04</b> Nylon fabric Colour: blue <b>TS04C B/R</b> Cordura fabric Colour: blue, red Wipe-off foil <b>TS04V B/R/G/Y</b> Colour: blue, red, green, yellow Length-Width-Height: 40-20-10cm Weight: 0,25kg</p>
<p>Carry bag for splint set</p> 	<p><b>TS06</b> Nylon fabric Colour: blue <b>TS06C B/R</b> Cordura fabric Colour: blue, red Wipe-off foil <b>TS06V B/R/G/Y</b> Colour: blue, red, green, yellow Length-Width-Height: 65-45-12cm Weight: 0,4kg</p>
<p>Carry bag for mattress with backpack straps</p> 	<p><b>TM01</b> Nylon fabric Colour: blue <b>TM01C B/R</b> Cordura fabric Colour: blue, red <b>TM01V B/R/G/Y</b> Wipe-off foil Colour: blue, red, green, yellow for mattress with width less than 90cm Length-Width-Height: 80-45-25cm Weight: 0,7kg</p>
	<p><b>TM02</b> Nylon fabric Colour: blue <b>TM02C B/R</b> Cordura fabric Colour: blue, red <b>TM02V B/R/G/Y</b> Wipe-off foil Colour: blue, red, green, yellow for mattress with width of more than 90cm Length-Width-Height: 100-45-30cm Weight: 1kg</p>

## Accessories for Vacuum Products

<p>Property indicator label</p>		<p>Label sizes: 8cm x 5cm 17cm x 7cm A4</p>
<p>9022 Repair Kit with self-adhesive patches</p>		



# TRS Thermo Rescue Sack

for use with a vacuum mattress  
and/or spine boards or rescue stretchers



- Water repellent and rugged material provides good protection without additional blanket
- Designed for easy access to the patient
- 10 carry handles
- Under surface with antiskid stripes
- Zipper and Hook-and-Loop closing for fast handling and easy access to the patient
- Optional property indicator label



Undersurface with antiskid stripes

Easy access to the patient

Fits also large Patient

### Specification

<b>TRSX01</b>	Material: Nylon, Carry Handles: 10 Width Head/Middle/Foot: 45cm/70cm/40cm Length: 200cm Weight: 2,5kg
<b>TRSX02</b>	Material: Nylon, Carry Handles: 10 Width Head/Middle/Foot: 45cm/70cm/40cm Length: 220cm Weight: 2,9kg
<b>TRSX03</b>	Material: CORDURA quilted, Carry Handles: 10 Width Head/Middle/Foot: 45cm/70cm/40cm Length: 220cm Weight: 2,9kg

# HRS Rescue Sack

for use with a vacuum mattress  
and/or spine boards or rescue stretchers



- Made of water repellent rugged CORDURA (HRS01) or water repellent heat insulating SOFT SHELL (HRS02)
- Designed for easy access to the patient
- 10 carry handles
- Under surface with antiskid stripes
- Zipper and Hook-and-Loop closing for fast handling and easy access to the patient
- Optional property indicator label
- CE marked



<b>HRS01</b>	Material: Cordura, red, Carry Handles: 10 Width Head/Middle/Foot: 60cm/60cm/40cm Length: 215cm Weight: 2,5kg
<b>HRS02</b>	Material: Softshell, orange, Carry Handles: 10 Width Head/Middle/Foot: 60cm/60cm/40cm Length: 215cm Weight: 2,5kg

# BRS Helicopter Rescue Sack

for use with a vacuum mattress  
and/or spine boards or rescue stretchers



- Optimised for rescue with helicopter and use with vacuum mattress
- The vacuum mattress can be deflated without opening the rescue bag
- Antiskid under surface to prevent slipping on snow
- Hook and loop fastener and zipper for easy access to the patient and to avoid wind-sack phenomenon and rotation
- Stable head protection
- Lifting strap system is adjustable to achieve horizontal position of the patient
- Patient restraint straps with metal buckles
- Pocket for vacuum pump
- Plated-through holes for free drainage of water
- Anti-rotation rings
- CE labelled and certified for helicopter lifting

## BRS03

All straps are sewn onto the rescue bag and are stowed in easy to open pockets



## BRS04

The detachable lifting bridle is attached to the carry handles with karabiner



Easy access to the patient



Integrated bag for the vacuum pump



Undersurface with antiskid stripes



## Specification



### BRS03

### BRS04

Carry Handles: 10  
Patient restraint straps: 5  
Width, Head-Middle-Foot: 40-70-45cm  
Length: 220 cm  
Weight: 4 kg

### TBRS03

Carry bag  
Width-Length-Height: 44-80-33 cm  
Weight: 1 kg

# KOHLBRAT & BUNZ

*Emergency & Rescue Systems*

*Made in Austria*

Quality management in accordance with  
EN ISO 9001 - EN ISO 13485

certified by



*made by*

**KOHLBRAT & BUNZ GMBH**



**A-5550 Radstadt**

**Austria**

**T: +43 6452 7193 0**

**F: +43 6452 7193 51**

**office@kohlbrat-bunz.com**

**www.kohlbrat-bunz.com**

**www.redvac.com**

Pictures are symbolic and may differ from the original. Stated weights are approximate.  
Subject to technical changes.